

cooltech

technical specifications

basic operation

type of equipment controlled cooling device for medical use

technical characteristics

minimum temperature - 8°C
handpieces Straight **HP**, Curved **HP**, Tight **HP**
Double **HP**, Tiny **HP**, Oval **HP**

power supply

power supply Red monofase, 100-240 V
maximum power 1000 W
cooling capacity in each handpiece 90 W programming 0 °C
with a Tamb of 25 °C

regulation

class of equipment as standard IEC 60601:2005 Class I, type BF
clasification according to directive 2007/47/CE Class II a
marking CE0051
the device complies with the requeriments Directive 2007/47 CE (93/42/CE)
of directives (Medical Device)
Directive 2002/95 CEE (RoHS)
Directive 2002/96 ECC (RAEE)
own patent EP-12382287

physical characteristics

dimensions 540 x 440 x 550 mm
weight 39 kg
screen touch screen 10,4"
