

# cooltech

## technical specifications

---

### basic operation

type of equipment controlled cooling device for medical use

---

### technical characteristics

minimum temperature - 8°C  
handpieces Straight **HP**, Curved **HP**, Tight **HP**  
Double **HP**, Tiny **HP**, Oval **HP**

### power supply

power supply Red monofase, 100-240 V  
maximum power 1000 W  
cooling capacity in each handpiece 90 W programming 0 °C  
with a Tamb of 25 °C

---

### regulation

class of equipment as standard IEC 60601:2005 Class I, type BF  
clasification according to directive 2007/47/CE Class II a  
marking CE0051  
the device complies with the requeriments Directive 2007/47 CE (93/42/CE)  
of directives (Medical Device)  
Directive 2002/95 CEE (RoHS)  
Directive 2002/96 ECC (RAEE)  
own patent EP-12382287

---

### physical characteristics

dimensions 540 x 440 x 550 mm  
weight 39 kg  
screen touch screen 10,4"

---